

schülke -†

A strong heart for CSSD circulation

Efficiency and cost-effectiveness
with the schülke pumping station



the plus of pure
performance

The schülke pumping station

compact – efficient – cost effective



For the optimal efficiency of central sterilisation it is becoming more common for the supplies for the cleaning and disinfection machines to be stored at a different location to the necessary cleaning chemicals.

The advantages are obvious: the process chemicals are supplied in large containers which do not need to be changed as often and the supplying of the machines does not disrupt the work processes in the CSSD (Central Sterile Supplies Department). In this way also any possible hazard to the staff resulting from handling the chemicals is minimised.

Many dosing systems require substantial space and offer complicated functions which are unnecessarily in daily work. schülke has developed the perfect solution for the CSSD.

The schülke pumping station offers robust technology “Made in Germany” in the smallest amount of space. Modern components can be found in a space-saving, quality high-grade steel casing which work with a minimum of maintenance and without any complicated programming.

As well as providing the chemicals schülke also offers a perfectly customised pumping station which supplies the installed automatic cleaning and disinfecting machines with the chemicals required for the process. The system pumps up to 6 different products and supplies the washing machines in the CSSD.

A continuous supply for the machines is ensured by intelligent, unmanned control mechanisms – without any programming input by the user.

The pumped amount, level regulating and function signals are preset by our technicians; the user can operate the system immediately and without any laborious training.



Your advantages:

The schülke pumping station is available in two basic versions (model 400 and model 600) and offers you the following features:



- The clearly laid out display panel of the control cabinet displays the status of the system and the individual pumps.
- Free of complicated control and operating elements.
- Flexible for possible expansion of the CSSD.
- Model 400 – 4 pumps | Model 600 – 6 pumps



- The display module quickly and clearly provides information on the operating status and empty-status or, respectively, full-status signals as well as possible faults.
- Optionally also available with an acoustic signal.



- The pump unit contains powerful pumps, which only need very little space due to their efficiency.
- Reliable technology, able to be operated with the minimum of maintenance requirements.
- Space-saving due to the highly-efficient feed pumps.

The chemical-resistant suction lances are equipped with an electronic level regulator; the need to replace the container is displayed in good time.



The reservoir container-tank is equipped with a fully electronic level regulator which ensures a continuous supply to the washing equipment.

Service from one source: your contact for further information and documents:

Customer service | Sales | Service

Phone 040/52100-214 Phone 040/52100-269
Fax 040/52100-228 Fax 040/52100-228

Postal address: 22840 Norderstedt

- 1 | **Dirk Bartels**
Email: dirk.bartels@schuelke.com
Mobile: 0171/2464119
- 2 | **Dirk Winter**
Email: dirk.winter@schuelke.com
Mobile: 0171/6379667
- 3 | **Tobias Kring**
Email: tobias.kring@schuelke.com
Mobile: 0172/9982449
- 4 | **Markus Radon**
Email: markus.radon@schuelke.com
Mobile: 0172/9264850
- 5 | **Heiner Franke**
Email: heiner.franke@schuelke.com
Mobile: 0171/2494003
- 6 | **Felix Grohs-Vekeman**
Email: felix.grohs-vekeman@schuelke.com
Mobile: 0171/6305470
- 7 | **Detlef Göddel**
Email: detlef.goeddel@schuelke.com
Mobile: 0172/4058022
- 8 | **Maik Böttcher**
Email: maik.boettcher@schuelke.com
Mobile: 0171/2466104
- 9 | **Wenzel Holecek**
Email: wenzel.holecek@schuelke.com
Mobile: 0171/2495594
- 10 | **Rene Schmidt**
Email: rene.schmidt@schuelke.com
Mobile: 0172/4058020
- 11 | **Stefan Kletzenbauer**
Email: stefan.kletzenbauer@schuelke.com
Mobile: 0171/6375820



Schülke & Mayr GmbH
22840 Norderstedt | Germany
Phone | Fax +49-40 52100-0 | -318
www.schuelke.com

A company of the
Air Liquide Group 